

## INSTRUCTIONS



TRIM SIZE:  
3.875" w X 9.5" h

### SANDAL CONVERSION KIT

Enjoy this sandal conversion kit as our free gift to you for purchasing WalkFit® Platinum orthotics.

## I M P O R T A N T

**READ ALL INSTRUCTION MATERIALS BEFORE USE OF THE WALKFIT® PLATINUM ORTHOTICS AND THIS SANDAL CONVERSION KIT. IF YOU HAVE ANY CONCERNS ABOUT WHETHER OR NOT YOU SHOULD USE YOUR ORTHOTICS WITH SANDALS, PLEASE CONSULT WITH YOUR PHYSICIAN.**


Now you can get the benefits of our orthotics in your favorite pair of sandals. The WalkFit® Platinum orthotics sandal conversion kit contains eight white circular pairs of hook and loop adhesives. Some of the hook and loop adhesive circles feel soft to the touch and some feel rough and prickly like sand paper. After you have located all the hook and loop adhesive circle pairs and separated the "rough" ones from the "soft" ones, proceed as follows:

**1** Peel the "rough" hook-side adhesive circles off of the strip that they came on and attach adhesive-side down to the bottoms of your orthotics – two near the front portion, and two near the back (as shown in the diagram below). You can always pull them off and re-adjust their positioning to find what works best for you.

**2** Next, place the "soft" loop-side of adhesive circles inside your sandals, coordinating their placement to line up with the "rough" circles attached to the bottoms of the WalkFit® Platinum orthotic. Make sure that all hook and loop circles are lined up and firmly attached before use to ensure your safety, comfort and optimal WalkFit® Platinum orthotic performance.

\*When your hook and loop adhesives start to wear off, replace immediately with new adhesive circles (usually available at craft and/or sewing stores nationwide).

#### DIAGRAM:



© 2004-2008 Walkfit, LLC. All rights reserved. U.S. Patent Nos. D509,951; 6,990,756. Chinese Patent No. 200430011865.5. European Registered Design No. 313317-0001. Hong Kong Patent No. 0412508.6. Additional patents pending. Made in China.

WFSC050908